

Sinusitis
Figure 1

FOREF 0251001






FIGURE 2

① Test kit




② Open the seal




③ Lay over nose

hole sealing fold

vent holes




④ Collect discharge




⑤ Spread with fingers

Reagent separators



⑥ Slide stick over discharge



⑦ Reading

	A	B	C	D	E
1					
2					
3					
4					
5					

SCORES

Score and chart

- A Leukocyte esterase
- B Nitrite
- C pH
- D Protein
- E TAME esterase

FIGURE 3